

# Enhancing monoclonal antibody production with picodroplet technology.

Tara-Jane Bloomfield-Gerdes, Giuseppe Benazzi, Karine Enesa  
Sphere Fluidics Limited, The Jonas Webb Building, Babraham Research Campus, Cambridge, CB22 3AT, UK

## Introduction

**Single cell cloning** is a crucial step in **monoclonal antibody production**. It can be achieved through costly, challenging and time consuming processes, e.g. FACS or limiting dilution.

We present a novel platform, **Cyto-Mine®**, that enables **ultra-high throughput sorting** of single cells and the **detection of their secreted products** in miniaturised (pL to nL) volumes called **picodroplets**. Using Chinese Hamster Ovary (CHO), we demonstrate that **Cyto-Mine®**:



1. Improves viability of single cells when encapsulated in picodroplets compared to traditional FACS sorting.
2. Allows accurate measurement of antibody secretion levels from encapsulated single cells within 1-2 hours.

There are many applications for this technology, including:

- **Biopharmaceutical discovery:** antibody (transcript) discovery from primary plasma cells, B-cells or hybridomas.
- **Bioprocessing:** enables measurement and isolation of clones that are high expressors of antibodies.

## Materials and methods


Flow diagram showing the encapsulation and sorting of a single cell:



To generate picodroplets containing a single cell, **Pico-Surf™ 1, 5% in Novec 7500** (surfactant dissolved in fluorinated oil) and cell suspension are injected into the channels of the picodroplet generation biochip. The flow rates are controlled by syringe pumps. To generate **300 pL** picodroplets, the flow rates for the cell suspension and oil are  $1000 \mu\text{L h}^{-1}$  and  $1400 \mu\text{L h}^{-1}$  respectively. With those conditions, we can generate **3.6 million picodroplets** per hour. The picodroplets are then sorted based on fluorescence, which correlates with the amount of antibody produced by each single cell.


## Results

**Figure 1: Poisson distribution of encapsulated CHO cells over time.**




**Figure 1:** Poisson distribution of CHO cells in 300 pL picodroplets over time compared to the theory. Cells were stained with  $0.5 \mu\text{M}$  Calcein AM, and 300 pL picodroplets were generated with aliquots collected at dedicated time points. Images were taken with a fluorescent microscope and analysed visually.

**Figure 2: Viability and recovery of CHO cells after sorting.**



**Figure 2:** Viability and recovery of CHO cells compared to FACS sorting. **A)** CHO cells were co-stained with  $0.2 \text{ mM}$  Cell Tracker™ Green CMFDA and  $3 \text{ mM}$  DRAQ7™. Stained cells underwent single cell sorting (as described in materials and methods), or were sorted using a BD FACSAria™ III cell sorter. Cell viability was assessed by detecting the ratio between green and red cells using fluorescent microscopy. **B)** Cells were not stained, and underwent single cell sorting (as described in materials and methods), or were sorted (and dispensed) using a BD FACSAria™ III cell sorter. The cells were then serially diluted to statistically give 1 cell per well, and plated on a 96 well plate. After approximately 4 weeks, visible clones were counted.

**Figure 3: Detection of antibody secretion from CHO cells in picodroplets.**



**Figure 3:** Long term viability of CHO cells in 300 pL picodroplets, and the amount of antibody secretion from CHO cells in 300 pL picodroplets compared to bulk. CHO cells were co-stained with  $0.2 \text{ mM}$  Cell Tracker™ Green CMFDA and  $3 \text{ mM}$  DRAQ7™. **A)** After picodroplet generation the encapsulated cells were incubated at  $37^\circ\text{C}$  for the indicated amount of time before viability was assessed. Cell viability was assessed by detecting the ratio between green and red cells using fluorescent microscopy. **B)** CHO cells were encapsulated in the presence of Goat anti-human IgG Fc DyLight 488 and Goat anti human (Fab')<sub>2</sub> DyLight 594. After picodroplet generation the encapsulated cells were incubated at  $37^\circ\text{C}$  for the indicated amount of time. Antibody secretion was measured using the SpectraMax® i3 Multi-Mode Detection Platform from Molecular Devices.

## Conclusion

We showed that Cyto-Mine® is an alternative technology to perform single cell cloning with statistically significantly higher recovery rates than flow cytometry. It can also be used for high-throughput screening of antibody-secreting cells enabling dramatic savings in consumable costs and improved assay sensitivity.

## References

- Chokkalingam *et al.*, 2013, *Lab Chip*, **13**, 4740-4744.  
Smith *et al.*, 2013, *Anal. Chem.*, **85** (8), 3812-3816.  
Mazutis *et al.*, 2013, *Nature Protocols*, **8**, 870-891.  
Debs *et al.*, 2012, *PNAS*, **109** (29), 11570-11575.

## Contact

Dr Karine Enesa (Cell & Molecular Biology):  
[Karine.Enesa@spherefluidics.com](mailto:Karine.Enesa@spherefluidics.com)

[www.spherefluidics.com](http://www.spherefluidics.com)  
Tel: +44 (0)1223 804200

