


ESI-Mine™

A label-free platform for high-throughput miniaturized electrospray ionization mass spectrometry (ESI-MS)



Recover rare cells producing proteins or other molecules of interest based on the result of the MS analysis

ESI-Mine™

ESI-Mine™ is a high-throughput, label-free platform for MS-based screening. It uniquely allows live clone retrieval of rare cells producing proteins or other molecules of interest based on the result of concurrent MS analysis.

Using our proprietary picodroplet (picolitre volume droplet) technology, ESI-Mine™ can test up to 200,000 biomolecular samples (e.g. peptides, enzymes, antibodies, metabolites and small molecules) per day in miniaturized volumes of 500-700pL (picodroplets) facilitating more efficient screening and analysis in bioproduction or synthetic biology processes.

Applying this novel approach, researchers can screen a large library of picodroplets containing microorganisms producing proteins or small molecules of interest. One major drawback of using MS in analysis is that samples are usually destroyed and cannot be retrieved. To overcome this total sample destruction by ESI-MS, we have developed a unique microfluidic workflow where picodroplets are split into two replicate streams with identical clones, one stream is analyzed by ESI-MS, and the second stream is subsequently 'held' and only those 'hit' daughter picodroplets sorted and retrieved based on the result of the MS analysis.



ESI-Mine™
Platform

Applications

ESI-Mine™ is suitable for several applications including:

Synthetic Biology: Facilitates study of valuable molecules produced by libraries of engineered microbes.

Bioprocessing: Enables analysis and identification of microbial clones that are high expressors of, e.g. Active Pharmaceutical Ingredients.

Drug and Molecular Interaction Studies: Accurately measures drug or drug fragment interaction with soluble proteins for improved structure-activity relationship (SAR) studies.

Enzyme Evolution/Site-Directed Mutagenesis: Helps determine the efficiency of engineered enzymes on their native substrates.


The Technology

MS offers label-free detection of the presence, absence, and abundance of multiple, specific molecules with high sensitivity and selectivity. However, it is relatively low-throughput and sample destruction during the analysis prevents the direct sorting of 'hit' molecules or samples of interest.

The microfluidic-based workflow of ESI-Mine™ addresses the main challenges in ESI-MS while offering advantages in high speed, increased sample throughput and lower cost over existing solutions.

Cells are encapsulated into picodroplets in a one-droplet-one-cell (OCOD) manner, creating a library of picodroplets that can express and retain novel proteins and/or other hydrophilic entities. The picodroplets are incubated to enable cell growth, protein expression and metabolic processes to occur and are then reinjected into the proprietary microfluidic emitter biochip on the ESI-Mine™ platform for label-free screening followed by real-time sorting based on MS analysis results. At this stage, each picodroplet is asymmetrically split, producing two daughter picodroplets from the same parent picodroplet. The larger daughter picodroplet is analyzed by ESI-MS. Its smaller sibling picodroplet is processed through microfluidic channels and sorted based on the MS signal. Hits are then collected as a pool for further downstream analysis.

Complete workflow of single cell analysis and sorting on ESI-Mine™ platform



Benefits

Recover Rare Variants

ESI-Mine™ uniquely allows live clone retrieval of rare cells producing proteins or other molecules of interest based on the result of the MS analysis.

High Throughput

ESI-Mine™ can selectively sort up to 200,000 biomolecular samples per day for efficient and high-throughput analysis of large cell libraries.

Miniaturization


Using picodroplets, the amount of media required is significantly reduced, further reducing costs.

Flexibility

We offer a range of customizable microfluidic biochips enabling the development of complex workflows on a single microfluidic device and interfacing to different MS instrument models

Versatility

ESI-Mine™ enables a broad range of novel microfluidic workflows and applications.



ESI-Mine™
Biochip

High-throughput and cost-effective bioprocessing


The table below summarises the substantial cost savings that can be achieved with the implementation of picodroplet technologies compared to a microplate format approach.


| Parameter | Microplate | Picodroplet |
|---------------------------------------|--|---|
| Throughput samples/day | 10 ² - 10 ³ (LC-MS) 10 ⁴ (RF-MS) | 10 ⁵ - 10 ⁶ |
| Reagent volume (to screen 1M mutants) | 100 L (e.g. £500,000 for 10 mM NADPH) | 1 mL (e.g. £5 for 10 mM NADPH) |
| Equipment | <ul style="list-style-type: none">• Qpix colony picker• Multi-plate Incubators• Centrifuges• FX/Tecan liquid handling robots• LC-MS or RF-MS• MTPs (£20,000-£50,000) | <ul style="list-style-type: none">• Syringe or air pressure pumps• Fabricated chips• Mass Spec• Biochip (£500). |

Specifications

| Run Specifications | |
|-----------------------------|--|
| Sample input method | Loaded into single-use disposable Picodroplet Reinjection & MS Emitter microfluidic biochip |
| Sample input format | Picodroplet |
| Workflows (operation modes) | Synthetic Biology; Enzyme Evolution; Drug Screening |
| Containment and sterility | Requires biological safety lab |
| Detection | |
| Detection system | TOF MS and Laser-Induced Fluorescence (e.g. fluorophores) |
| Excitation wavelengths | 488 nm |
| Fluorescent wavelengths * | Below 506nm and 524-634nm |
| Camera | High-speed CMOS |
| PC | |
| Computer | Embedded internally as part of the TOF MS; Control PC |
| PC operating system | Microsoft Windows 10 Professional |
| Monitor | Colour LCD (24") |
| Software Specifications | |
| System control software | TOF MS Driver; ESI-Mine™ Studio; Elveflow® ESI Pressure control software; StreamPix 6; Coherent Connection |
| Data tracking | On screen data point |
| Work Environment | |
| Clearance | 30 cm |
| Humidity | 30 - 80% |
| Operating temperature | 16°C - 26°C |
| Site preparation | See the ESI-Mine™ System Site Requirements Guide |
| Consumables | |
| Microfluidic biochips | Picodroplet Reinjection & MS Emitter microfluidic biochip; Picodroplet Reinjection, Splitting, MS Emitter & Sorting microfluidic biochip |
| Specialist chemicals | Picodroplet Oil; Picodroplet Reinjection Oil; Picodroplet Spacing Oil |

Implementation Example: ESI-Mine™ Linked to a Perkin Elmer Axion 2





ESI-Mine™

A label-free platform for high-throughput miniaturized electrospray ionization mass spectrometry (ESI-MS)

| Code | Product Ordering Information |
|------|------------------------------|
| S004 | ESI-Mine™ |

Notes:

Sphere Fluidics Ltd is an ISO 9001:2015 accredited company. Pico-Surf™, Pico-Break™, Pico-Glide™, Pico-Wave™, Pico-Gen™, Pico-Sort™, Pico-Mix™, Pico-Safe™ are trademarks and/or patented technologies of Sphere Fluidics Ltd.

All Sphere Fluidics' supplied chemicals and bioreagents are Animal Origin Free and GLP-compliant.

For research and development purposes only.

Product specifications subject to change without notice.

Content ©Sphere Fluidics Ltd.



Contact us

Email: Info@spherefluidics.com

Service: Service@spherefluidics.com

Web: www.spherefluidics.com

Visit our website for a current list of our distributors:

www.spherefluidics.com/contact-us/distributors/